

Columbia Coatings

Technical Data Sheet

Columbia Coatings
1173 Industrial Park Rd
Columbia, TN 38401
931-388-7730
931-388-5573 FAX



Number:	CPA-BCCS
Name:	Big Chimpin' Chemical Stripper
Type:	Liquid Stripper

Recommended Application

Step #1	Prepare a strip tank with lid and proper ventilation.
Step #2	Submerge parts in heated solution
Step #3	Recommended dwell time is 30-60 minutes in heated bath (dwell time may vary). For best results heat solution to 175-185 °F
Step #4	Rinse parts using fresh water and pressure washer.

Benefits

- Especially designed for magnesium, aluminum and non-ferrous metals. Can be used on ferrous metals.
- No methylene chloride- safe for use.
- Low vapor pressure reduces VOC's, Environmentally friendly

Health and Safety

- Please refer to SDS for complete health and safety information.

Storage

- Store in a dry location at room temperature. Keep container closed and maintain all original markings and labels. Keep this container and vapors from this container away from heat, sparks, flame, and other ignition sources. Recommended storage in ferrous, stainless steel and plastic containers only.

www.columbiacoatings.com